



## ClimbUP® (the original)

The ClimbUP® makes bed bug monitoring easy. When used under furniture legs where people are sleeping, the ClimbUP® will make visible very low levels of bed bug infestation to any observer. The exterior climbing surface is white paper. The ClimbUP® is round with total diameter of 6" and height of 1". The center well diameter is 3 3/4".



## ClimbUP® BG (Black Grip)

The ClimbUP® BG is used in the same manner as Original ClimbUP®. The exterior climbing surface is black paper. The ClimbUP® BG is round with total diameter of 6" and height of 1". The center well diameter is 3 3/4".



## ClimbUP® HD (Hotel Discreet)

The ClimbUP®HD is small, black and discreet. Used under standard metal bed frame casters and glides, it appears to be a coaster and it hides captured bed bugs until swabbed out. The exterior climbing surface of the Climbup®HD is black paper. The Climbup®HD is round with total diameter of 4" and height of 1". The center well diameter is 2 1/2".



## ClimbUP® XL

The ClimbUP®XL can be used under larger furniture legs or double headboard legs to detect bed bug infestations. The exterior climbing surface is white paper. The Climbup®XL is an 8" x 10" rectangle with a height of 2". The center well is 4" x 6".



## ClimbUP® T

The ClimbUP®T can be used in an unoccupied room with a CO<sub>2</sub> Tank or insulated beverage cooler of dry ice in the center well or in an occupied room under very large furniture legs to detect bed bugs. The exterior climbing surface of Climbup®T is black paper. The Climbup®T is round with total diameter of 9 1/4" and height of 2". The center well diameter is 5 3/4".